



SPINE SOLUTION



ORTHOWIND

SPINE SOLUTION



about us

Orthowind, our mission is to provide high-quality spine, trauma and surgical instruments that enhance patient care and improve surgical outcomes. We are dedicated to supporting healthcare professionals with precision-engineered tools that meet the highest standards of safety, innovation, and reliability. By prioritizing excellence in both product and service, we strive to ensure that surgeons and medical teams can deliver the best possible.

Orthowind Medical

phone +90 541 835 5589

address Sarıbeyler Mah. Sarıbeyler Küme Evler No:
462 Keçiören / Ankara / Türkiye

email info@orthowind.com.tr

website www.orthowind.com.tr



Contents

Cervical Spine Products 1.0 - 1.15

Cervical Bladed Peek Cage	1.2
Cervical Peek Cage	1.4
Expandable Cervical Peek Cage	1.6
Mobile Cervical Disc Prosthesis	1.8
Cervical Plate and Screws	1.10
Cervical Corpectomy Cage	1.12
Cervical Spinal System	1.14

Lumbar Spine Products 2.0 - 2.22


Spinal System 6.0	2.2
Polar Spinal System 5.5 / 6.0	2.8
Lumbar Corpectomy Cage	2.19
Interspinous Device	2.20
Tlif Peek Cage	2.22
Plif Peek Cage	2.24
Expandable Plif Peek Cage	2.26

Pedis Spinal System 3.0 - 3.4



Cervical Spine Products

Cervical Bladed Peek Cage



Features

- Cervical Bladed Peek Cage is manufactured by using PEEK and Titanium, which is compatible with MRI and CT and which does not result in permanent lesions.
- Blades for a more reliable holding between the endplates.
- Anterior cervical plate may not required for supplemental fixation.





Code	Height	Lenght	Width
ORT-MCPCB41214	4	12	14
ORT-MCPCB51214	5	12	14
ORT-MCPCB61214	6	12	14
ORT-MCPCB71214	7	12	14
ORT-MCPCB81214	8	12	14


Code	Height	Lenght	Width
ORT-MCPCB41216	4	12	16
ORT-MCPCB51216	5	12	16
ORT-MCPCB61216	6	12	16
ORT-MCPCB71216	7	12	16
ORT-MCPCB81216	8	12	16

Code	Height	Lenght	Width
ORT-MCPCB41414	4	14	14
ORT-MCPCB51414	5	14	14
ORT-MCPCB61414	6	14	14
ORT-MCPCB71414	7	14	14
ORT-MCPCB81414	8	14	14

Code	Height	Lenght	Width
ORT-MCPCB41416	4	14	16
ORT-MCPCB51416	5	14	16
ORT-MCPCB61416	6	14	16
ORT-MCPCB71416	7	14	16
ORT-MCPCB81416	8	14	16




Cervical Peek Cage



Features

- Cervical Peek is manufactured by using PEEK material, which is compatible with MRI and CT and which does not result in permanent lesions.
- Easy insert and fit anatomic structure.
- Tantalum marker.
- Strong fixation with threaded surface and two titanium pins.






Code	Height	Lenght	Width
ORT-CPC041214	4	12	14
ORT-CPC051214	5	12	14
ORT-CPC061214	6	12	14
ORT-CPC071214	7	12	14
ORT-CPC081214	8	12	14

Code	Height	Lenght	Width
ORT-CPC041216	4	12	16
ORT-CPC051216	5	12	16
ORT-CPC061216	6	12	16
ORT-CPC071216	7	12	16
ORT-CPC081216	8	12	16

Code	Height	Lenght	Width
ORT-CPC041414	4	14	14
ORT-CPC051414	5	14	14
ORT-CPC061414	6	14	14
ORT-CPC071414	7	14	14
ORT-CPC081414	8	14	14

Code	Height	Lenght	Width
ORT-CPC041416	4	14	16
ORT-CPC051416	5	14	16
ORT-CPC061416	6	14	16
ORT-CPC071416	7	14	16
ORT-CPC081416	8	14	16




Expandable Cervical Peek Cage




Features

- Expandable Cervical Peek is manufactured by using PEEK material and expanding mechanism: Ti6Al4V- Eli alloy which is compatible with MRI and CT and which does not result in permanent lesions.
- Titanium pins together with titanium expanding mechanism let us identify the actual cage position.





Code	Height	Lenght	Width
ORT-ECPC125	5	12	14
ORT-ECPC126	6	12	14
ORT-ECPC127	7	12	14
ORT-ECPC128	8	12	14




Mobile Cervical Disc Prosthesis



Features

- Manufactured from Ti6Al4V-ELI material, which is compatible with MRI and does not allow any lesional problems.
- Superior and inferior surfaces of a toothed structure provides retention, so the prosthesis can be placed firmly in the upper and lower parts of the end plates.
- The surface provides safety with the fixation into superior and inferior vertebrae. ($\pm$) 10° extension/flexion and lateral bending, and ($\pm$) 8° rotation of the disc can be fixed/provided by a single tool.





Code	Height	x	Lenght	x	Width
ORT-CD12-04	4	x	12	x	14
ORT-CD12-05	5	x	12	x	14
ORT-CD12-06	6	x	12	x	14
ORT-CD12-07	7	x	12	x	14
ORT-CD12-08	8	x	12	x	14
ORT-CD12-09	9	x	12	x	14

Code	Height	x	Lenght	x	Width
ORT-CD14-04	4	x	14	x	14
ORT-CD14-05	5	x	14	x	14
ORT-CD14-06	6	x	14	x	14
ORT-CD14-07	7	x	14	x	14
ORT-CD14-08	8	x	14	x	14
ORT-CD14-09	9	x	14	x	14

Code	Height	x	Lenght	x	Width
ORT-CD16-05	5	x	14	x	16
ORT-CD16-06	6	x	14	x	16
ORT-CD16-07	7	x	14	x	16
ORT-CD16-08	8	x	14	x	16
ORT-CD16-09	9	x	14	x	16




L



H

W

Cervical Anterior Plate and Screws



Features

- The locking system that secures screws for one-step locking.
- The plate has a low profile and smooth surface designed to help minimize irritation.
- Simple surgical technique.



Code	Size	Level
ORT-CAP17	17 mm	1
ORT-CAP20	20 mm	1
ORT-CAP23	23 mm	1
ORT-CAP25	25 mm	1
ORT-CAP27	27 mm	2
ORT-CAP30	30 mm	2
ORT-CAP33	33 mm	2
ORT-CAP36	36 mm	2
ORT-CAP40	40 mm	2
ORT-CAP45	45 mm	3
ORT-CAP50	50 mm	3
ORT-CAP55	55 mm	3
ORT-CAP60	60 mm	3
ORT-CAP65	65 mm	3
ORT-CAP70	70 mm	4
ORT-CAP75	75 mm	4
ORT-CAP80	80 mm	4
ORT-CAP85	85 mm	4
ORT-CAP90	90 mm	4



Code	Size
ORT-CAPSM12	Ø3.5x12 mm
ORT-CAPSM14	Ø3.5x14 mm
ORT-CAPSM16	Ø3.5x16 mm
ORT-CAPSM18	Ø3.5x18 mm
ORT-CAPSM20	Ø3.5x20 mm



Code	Size
ORT-CAPS12	Ø4.0x12 mm
ORT-CAPS14	Ø4.0x14 mm
ORT-CAPS16	Ø4.0x16 mm
ORT-CAPS18	Ø4.0x18 mm
ORT-CAPS20	Ø4.0x20 mm



Code	Size
ORT-CAPSR12	Ø4.5x12 mm
ORT-CAPSR14	Ø4.5x14 mm
ORT-CAPSR16	Ø4.5x16 mm
ORT-CAPSR18	Ø4.5x18 mm
ORT-CAPSR20	Ø4.5x20 mm



Cervical Corpectomy Cage



Features

- Full contact with angled surface.
- Teeth on the surface, minimizing the risk of expulsion.
- Angled inferior and superior area allow a complete contact with vertebral surface and composed by one piece.
- With an efficient grafting space, the system allows applying graft before distraction and provides a one stage locking mechanism.





Code	Dia	x	Closed Lenght	x	Open Lenght	Angled
ORTCTC1216256	12	x	16	x	25	6°
ORTCTC1224406	12	x	24	x	40	6°
ORTCTC1239656	12	x	39	x	65	6°

Code	Dia	x	Closed Lenght	x	Open Lenght	Angled
ORTCTC1216256	14	x	16	x	25	6°
ORTCTC1224406	14	x	24	x	40	6°
ORTCTC1239656	14	x	39	x	65	6°

Code	Dia	x	Closed Lenght	x	Open Lenght
ORTCTC101013	10	x	10	x	13
ORTCTC101317	10	x	13	x	17
ORTCTC101625	10	x	16	x	25

Code	Dia	x	Closed Lenght	x	Open Lenght
ORTCTC121013	12	x	10	x	13
ORTCTC121317	12	x	13	x	17
ORTCTC121625	12	x	16	x	25
ORTCTC122440	12	x	24	x	40
ORTCTC123965	12	x	39	x	65

Code	Dia	x	Closed Lenght	x	Open Lenght
ORTCTC141013	14	x	10	x	13
ORTCTC141317	14	x	13	x	17
ORTCTC141625	14	x	16	x	25
ORTCTC142440	14	x	24	x	40
ORTCTC143965	14	x	39	x	65

Code	Dia	x	Closed Lenght	x	Open Lenght
ORTCTC161013	16	x	10	x	13
ORTCTC161317	16	x	13	x	17
ORTCTC161625	16	x	16	x	25
ORTCTC162440	16	x	24	x	40
ORTCTC163965	16	x	39	x	65

Code	Dia	x	Closed Lenght	x	Open Lenght
ORTCTC201013	20	x	10	x	13
ORTCTC201317	20	x	13	x	17
ORTCTC201625	20	x	16	x	25




Cervical Corpectomy Cage Screw

Code
ORTCCLS

Cervical

Spinal System


Colour Parameters as per Ø



Turquoise
Ø 3.5



Pink
Ø 4.0

Golden
Ø 4.5



Features


- Self-tapping flute centers screw for easy insertion.
- Low profile screw head 3 diameters with cancellous profile (3.5 mm & 4.0 mm and 4.5 mm Rescue Screw) 3.5 mm diameter rod.
- Multiple screw, hook, rod, connector and plate options to better match patient anatomy.



Cervical Posterior Polyaxial Screw



Code	Size
ORT-PCPS3510	Ø3.5x10 mm
ORT-PCPS3512	Ø3.5x12 mm
ORT-PCPS3514	Ø3.5x14 mm
ORT-PCPS3516	Ø3.5x16 mm
ORT-PCPS3518	Ø3.5x18 mm
ORT-PCPS3520	Ø3.5x20 mm
ORT-PCPS3522	Ø3.5x22 mm
ORT-PCPS3524	Ø3.5x24 mm
ORT-PCPS3526	Ø3.5x26 mm
ORT-PCPS3528	Ø3.5x28 mm
ORT-PCPS3530	Ø3.5x30 mm
ORT-PCPS3532	Ø3.5x32 mm
ORT-PCPS4010	Ø4.0x10 mm
ORT-PCPS4012	Ø4.0x12 mm
ORT-PCPS4014	Ø4.0x14 mm
ORT-PCPS4016	Ø4.0x16 mm
ORT-PCPS4018	Ø4.0x18 mm
ORT-PCPS4020	Ø4.0x20 mm
ORT-PCPS4022	Ø4.0x22 mm
ORT-PCPS4024	Ø4.0x24 mm
ORT-PCPS4026	Ø4.0x26 mm
ORT-PCPS4028	Ø4.0x28 mm
ORT-PCPS4030	Ø4.0x30 mm
ORT-PCPS4032	Ø4.0x32 mm
ORT-PCPS4510	Ø4.5x10 mm
ORT-PCPS4512	Ø4.5x12 mm
ORT-PCPS4514	Ø4.5x14 mm
ORT-PCPS4516	Ø4.5x16 mm
ORT-PCPS4518	Ø4.5x18 mm
ORT-PCPS4520	Ø4.5x20 mm
ORT-PCPS4522	Ø4.5x22 mm
ORT-PCPS4524	Ø4.5x24 mm
ORT-PCPS4526	Ø4.5x26 mm
ORT-PCPS4528	Ø4.5x28 mm
ORT-PCPS4530	Ø4.5x30 mm
ORT-PCPS4532	Ø4.5x32 mm



Cervical Posterior Polyaxial Screw (Half Thread)

Code	Size
ORT-PCPSC3510	Ø3.5x10 mm
ORT-PCPSC3512	Ø3.5x12 mm
ORT-PCPSC3514	Ø3.5x14 mm
ORT-PCPSC3516	Ø3.5x16 mm
ORT-PCPSC3518	Ø3.5x18 mm
ORT-PCPSC3520	Ø3.5x20 mm

Cervical Rod 3.5

Code	Size
ORT-SR1304	Ø3.5x40 mm
ORT-SR1308	Ø3.5x80 mm
ORT-SR1312	Ø3.5x120 mm
ORT-SR1320	Ø3.5x200 mm
ORT-SR1324	Ø3.5x240 mm
ORT-SR1345	Ø3.5x(200-70) mm

Cervical Lateral Connector

Code	Size
ORT-LCNT12	10 mm
ORT-LCNT18	18 mm

Cervical Straight Transverse Connector

Code	Size
ORT-TLC130	30 mm
ORT-TLC135	35 mm
ORT-TLC140	40 mm
ORT-TLC145	45 mm
ORT-TLC150	50 mm
ORT-TLC155	55 mm



Cervical Transverse Link

Code	Size
ORT-TLRC130	30 mm
ORT-TLRC135	35 mm
ORT-TLRC140	40 mm
ORT-TLRC145	45 mm
ORT-TLHC	



Cervical Hook

Code	Size
ORT-PCH55	Ø4.5X5 mm
ORT-PCHM6	Ø4.5X6 mm
ORT-PCHL7	Ø4.5X7 mm



Cervical Lateral Connector

Code	Size
ORT-SDCN02	For 1 Rod
ORT-SECN01	For 2 Rod



Rigid Plate

Code	Size
ORT-CPOPL	L
ORT-CPOPM	M
ORT-CPOPS	S

Rigid Plate Screw

Code	Size
ORT-CPOSS4	Ø4.0x4 mm
ORT-CPOSS5	Ø4.0x5 mm
ORT-CPOSS6	Ø4.0x6 mm
ORT-CPOSL4	Ø4.5x4 mm
ORT-CPOSL5	Ø4.5x5 mm
ORT-CPOSL6	Ø4.5x6 mm



Set Screw

Code	
ORT-CSS	







Lumbar Spine Products

PROCYON

Spinal System 6.0mm Rod

Colour Parameters as per Ø

						
Pink Ø 4.5	Purple Ø 5.0	Blue Ø 5.5	Turquoise Ø 6.0	Golden Ø 6.5	Ice Blue Ø 7.0	Green Ø 7.5



Features

- Easy Lock System.
- Cylindrical feature Self-Tapping.
- 57 degree with polyaxial screws.
- Implantable pedicle screws as a Monoaxial, Polyaxial Cannulated screws.
- Pre-bent rods, rod types, hooks in different sizes, easy to use hand tools.
- Compatible with implants.
- More reliable tightening with the torx design of setscrew.
- Reverse angled setscrew thread design.

Polyaxial Screw



Code	Size
ORT-PAS3525	Ø3.5x25 mm
ORT-PAS3530	Ø3.5x30 mm
ORT-PAS3535	Ø3.5x35 mm
ORT-PAS3540	Ø3.5x40 mm
ORT-PAS3545	Ø3.5x45 mm
ORT-PAS4020	Ø4.0x20 mm
ORT-PAS4025	Ø4.0x25 mm
ORT-PAS4030	Ø4.0x30 mm
ORT-PAS4035	Ø4.0x35 mm
ORT-PAS4040	Ø4.0x40 mm
ORT-PAS4045	Ø4.0x45 mm
ORT-PAS4520	Ø4.5x20 mm
ORT-PAS4525	Ø4.5x25 mm
ORT-PAS4530	Ø4.5x30 mm
ORT-PAS4535	Ø4.5x35 mm
ORT-PAS4540	Ø4.5x40 mm
ORT-PAS4545	Ø4.5x45 mm
ORT-PAS5030	Ø5.0x30 mm
ORT-PAS5035	Ø5.0x35 mm
ORT-PAS5040	Ø5.0x40 mm
ORT-PAS5045	Ø5.0x45 mm
ORT-PAS5050	Ø5.0x50 mm
ORT-PAS5055	Ø5.0x55 mm
ORT-PAS5060	Ø5.0x60 mm
ORT-PAS5520	Ø5.5x20 mm
ORT-PAS5525	Ø5.5x25 mm
ORT-PAS5530	Ø5.5x30 mm
ORT-PAS5535	Ø5.5x35 mm
ORT-PAS5540	Ø5.5x40 mm
ORT-PAS5545	Ø5.5x45 mm
ORT-PAS5550	Ø5.5x50 mm
ORT-PAS5555	Ø5.5x55 mm
ORT-PAS5560	Ø5.5x60 mm

Code	Size
ORT-PAS6030	Ø6.0x30 mm
ORT-PAS6035	Ø6.0x35 mm
ORT-PAS6040	Ø6.0x40 mm
ORT-PAS6045	Ø6.0x45 mm
ORT-PAS6050	Ø6.0x50 mm
ORT-PAS6055	Ø6.0x55 mm
ORT-PAS6060	Ø6.0x60 mm
ORT-PAS6530	Ø6.5x30 mm
ORT-PAS6535	Ø6.5x35 mm
ORT-PAS6540	Ø6.5x40 mm
ORT-PAS6545	Ø6.5x45 mm
ORT-PAS6550	Ø6.5x50 mm
ORT-PAS6555	Ø6.5x55 mm
ORT-PAS6560	Ø6.5x60 mm
ORT-PAS7035	Ø7.0x35 mm
ORT-PAS7040	Ø7.0x40 mm
ORT-PAS7045	Ø7.0x45 mm
ORT-PAS7050	Ø7.0x50 mm
ORT-PAS7055	Ø7.0x55 mm
ORT-PAS7060	Ø7.0x60 mm
ORT-PAS7530	Ø7.5x30 mm
ORT-PAS7535	Ø7.5x35 mm
ORT-PAS7540	Ø7.5x40 mm
ORT-PAS7545	Ø7.5x45 mm
ORT-PAS7550	Ø7.5x50 mm
ORT-PAS7555	Ø7.5x55 mm
ORT-PAS7560	Ø7.5x60 mm
ORT-PAS8030	Ø8.0x30 mm
ORT-PAS8035	Ø8.0x35 mm
ORT-PAS8040	Ø8.0x40 mm
ORT-PAS8045	Ø8.0x45 mm
ORT-PAS8050	Ø8.0x50 mm
ORT-PAS8055	Ø8.0x55 mm

Code	Size
ORT-PAS8060	Ø8.0x60 mm
ORT-PAS8070	Ø8.0x70 mm
ORT-PAS8080	Ø8.0x80 mm
ORT-PAS8090	Ø8.0x90 mm
ORT-PAS80100	Ø8.0x100 mm
ORT-PAS8530	Ø8.5x30 mm
ORT-PAS8535	Ø8.5x35 mm
ORT-PAS8540	Ø8.5x40 mm
ORT-PAS8545	Ø8.5x45 mm
ORT-PAS8550	Ø8.5x50 mm
ORT-PAS8555	Ø8.5x55 mm
ORT-PAS8560	Ø8.5x60 mm
ORT-PAS8570	Ø8.5x70 mm
ORT-PAS8580	Ø8.5x80 mm
ORT-PAS8590	Ø8.5x90 mm
ORT-PAS85100	Ø8.5x100 mm






Code	Size
ORT-MAS3525	Ø3.5x25 mm
ORT-MAS3530	Ø3.5x30 mm
ORT-MAS3535	Ø3.5x35 mm
ORT-MAS3540	Ø3.5x40 mm
ORT-MAS3545	Ø3.5x45 mm
ORT-MAS4025	Ø4.0x25 mm
ORT-MAS4030	Ø4.0x30 mm
ORT-MAS4035	Ø4.0x35 mm
ORT-MAS4040	Ø4.0x40 mm
ORT-MAS4045	Ø4.0x45 mm
ORT-MAS4525	Ø4.5x25 mm
ORT-MAS4530	Ø4.5x30 mm
ORT-MAS4535	Ø4.5x35 mm
ORT-MAS4540	Ø4.5x40 mm
ORT-MAS4545	Ø4.5x45 mm
ORT-MAS5030	Ø5.0x30 mm
ORT-MAS5035	Ø5.0x35 mm
ORT-MAS5040	Ø5.0x40 mm
ORT-MAS5045	Ø5.0x45 mm
ORT-MAS5050	Ø5.0x50 mm
ORT-MAS5530	Ø5.5x30 mm
ORT-MAS5535	Ø5.5x35 mm
ORT-MAS5540	Ø5.5x40 mm
ORT-MAS5545	Ø5.5x45 mm
ORT-MAS5550	Ø5.5x50 mm
ORT-MAS5555	Ø5.5x55 mm

Code	Size
ORT-MAS6035	Ø6.0x35 mm
ORT-MAS6040	Ø6.0x40 mm
ORT-MAS6045	Ø6.0x45 mm
ORT-MAS6050	Ø6.0x50 mm
ORT-MAS6055	Ø6.0x55 mm
ORT-MAS6530	Ø6.5x30 mm
ORT-MAS6535	Ø6.5x35 mm
ORT-MAS6540	Ø6.5x40 mm
ORT-MAS6545	Ø6.5x45 mm
ORT-MAS6550	Ø6.5x50 mm
ORT-MAS6555	Ø6.5x55 mm
ORT-MAS7035	Ø7.0x35 mm
ORT-MAS7040	Ø7.0x40 mm
ORT-MAS7045	Ø7.0x45 mm
ORT-MAS7050	Ø7.0x50 mm
ORT-MAS7055	Ø7.0x55 mm
ORT-MAS7535	Ø7.5x35 mm
ORT-MAS7540	Ø7.5x40 mm
ORT-MAS7545	Ø7.5x45 mm
ORT-MAS7550	Ø7.5x50 mm
ORT-MAS7555	Ø7.5x55 mm
ORT-MAS8035	Ø8.0x35 mm
ORT-MAS8040	Ø8.0x40 mm
ORT-MAS8045	Ø8.0x45 mm
ORT-MAS8050	Ø8.0x50 mm
ORT-MAS8055	Ø8.0x55 mm

Polyaxial Spondylolisthesis Screw



Code	Size
ORT-PRS4530	Ø4.5x30 mm
ORT-PRS4535	Ø4.5x35 mm
ORT-PRS4540	Ø4.5x40 mm
ORT-PRS4545	Ø4.5x45 mm
ORT-PRS5530	Ø5.5x30 mm
ORT-PRS5535	Ø5.5x35 mm
ORT-PRS5540	Ø5.5x40 mm
ORT-PRS5545	Ø5.5x45 mm
ORT-PRS5550	Ø5.5x50 mm
ORT-PRS5555	Ø5.5x55 mm
ORT-PRS6035	Ø6.0x35 mm
ORT-PRS6040	Ø6.0x40 mm
ORT-PRS6045	Ø6.0x45 mm
ORT-PRS6050	Ø6.0x50 mm
ORT-PRS6530	Ø6.5x30 mm
ORT-PRS6535	Ø6.5x35 mm
ORT-PRS6540	Ø6.5x40 mm
ORT-PRS6545	Ø6.5x45 mm
ORT-PRS6550	Ø6.5x50 mm
ORT-PRS6555	Ø6.5x55 mm
ORT-PRS7035	Ø7.0x35 mm
ORT-PRS7040	Ø7.0x40 mm
ORT-PRS7045	Ø7.0x45 mm
ORT-PRS7050	Ø7.0x50 mm
ORT-PRS7055	Ø7.0x55 mm
ORT-PRS7535	Ø7.5x35 mm
ORT-PRS7540	Ø7.5x40 mm
ORT-PRS7545	Ø7.5x45 mm
ORT-PRS7550	Ø7.5x50 mm
ORT-PRS7555	Ø7.5x55 mm

Monoaxial Spondylolisthesis Screw



Code	Size
ORT-MRS4530	Ø4.5x30 mm
ORT-MRS4535	Ø4.5x35 mm
ORT-MRS4540	Ø4.5x40 mm
ORT-MRS4545	Ø4.5x45 mm
ORT-MRS5535	Ø5.5x35 mm
ORT-MRS5540	Ø5.5x40 mm
ORT-MRS5545	Ø5.5x45 mm
ORT-MRS5550	Ø5.5x50 mm
ORT-MRS6035	Ø6.0x35 mm
ORT-MRS6040	Ø6.0x40 mm
ORT-MRS6045	Ø6.0x45 mm
ORT-MRS6050	Ø6.0x50 mm
ORT-MRS6535	Ø6.5x35 mm
ORT-MRS6540	Ø6.5x40 mm
ORT-MRS6545	Ø6.5x45 mm
ORT-MRS6550	Ø6.5x50 mm
ORT-MRS7035	Ø7.0x35 mm
ORT-MRS7040	Ø7.0x40 mm
ORT-MRS7045	Ø7.0x45 mm
ORT-MRS7050	Ø7.0x50 mm
ORT-MRS7055	Ø7.0x55 mm
ORT-MRS7540	Ø7.5x40 mm
ORT-MRS7545	Ø7.5x45 mm
ORT-MRS7550	Ø7.5x50 mm
ORT-MRS7555	Ø7.5x55 mm



Code	Size
ORT-CPS4530	Ø4.5x30 mm
ORT-CPS4535	Ø4.5x35 mm
ORT-CPS4540	Ø4.5x40 mm
ORT-CPS4545	Ø4.5x45 mm
ORT-CPS4550	Ø4.5x50 mm
ORT-CPS4555	Ø4.5x55 mm
ORT-CPS5035	Ø5.0x35 mm
ORT-CPS5040	Ø5.0x40 mm
ORT-CPS5045	Ø5.0x45 mm
ORT-CPS5050	Ø5.0x50 mm
ORT-CPS5055	Ø5.0x55 mm
ORT-CPS5530	Ø5.5x30 mm
ORT-CPS5535	Ø5.5x35 mm
ORT-CPS5540	Ø5.5x40 mm
ORT-CPS5545	Ø5.5x45 mm
ORT-CPS5550	Ø5.5x50 mm
ORT-CPS5555	Ø5.5x55 mm
ORT-CPS6035	Ø6.0x35 mm
ORT-CPS6040	Ø6.0x40 mm
ORT-CPS6045	Ø6.0x45 mm
ORT-CPS6050	Ø6.0x50 mm
ORT-CPS6055	Ø6.0x55 mm
ORT-CPS6530	Ø6.5x30 mm

Code	Size
ORT-CPS6535	Ø6.5x35 mm
ORT-CPS6540	Ø6.5x40 mm
ORT-CPS6545	Ø6.5x45 mm
ORT-CPS6550	Ø6.5x50 mm
ORT-CPS6555	Ø6.5x55 mm
ORT-CPS7030	Ø7.0x30 mm
ORT-CPS7035	Ø7.0x35 mm
ORT-CPS7040	Ø7.0x40 mm
ORT-CPS7045	Ø7.0x45 mm
ORT-CPS7050	Ø7.0x50 mm
ORT-CPS7055	Ø7.0x55 mm
ORT-CPS7530	Ø7.5x30 mm
ORT-CPS7535	Ø7.5x35 mm
ORT-CPS7540	Ø7.5x40 mm
ORT-CPS7545	Ø7.5x45 mm
ORT-CPS7550	Ø7.5x50 mm
ORT-CPS7555	Ø7.5x55 mm
ORT-CPS8035	Ø8.0x35 mm
ORT-CPS8040	Ø8.0x40 mm
ORT-CPS8045	Ø8.0x45 mm
ORT-CPS8050	Ø8.0x50 mm
ORT-CPS8055	Ø8.0x55 mm

Monoaxial Cannulated And Cemented Screw

Code	Size
ORT-CMS5530	Ø5.5x30 mm
ORT-CMS5535	Ø5.5x35 mm
ORT-CMS5540	Ø5.5x40 mm
ORT-CMS5545	Ø5.5x45 mm
ORT-CMS5550	Ø5.5x50 mm
ORT-CMS5555	Ø5.5x55 mm
ORT-CMS6530	Ø6.5x30 mm
ORT-CMS6535	Ø6.5x35 mm
ORT-CMS6540	Ø6.5x40 mm
ORT-CMS6545	Ø6.5x45 mm
ORT-CMS6550	Ø6.5x50 mm
ORT-CMS6555	Ø6.5x55 mm
ORT-CMS7530	Ø7.5x30 mm
ORT-CMS7535	Ø7.5x35 mm
ORT-CMS7540	Ø7.5x40 mm
ORT-CMS7545	Ø7.5x45 mm
ORT-CMS7550	Ø7.5x50 mm
ORT-CMS7555	Ø7.5x55 mm



Rod 6.0




Code	Size
ORT-SR1604	Ø6.0x40 mm
ORT-SR1605	Ø6.0x50 mm
ORT-SR1606	Ø6.0x60 mm
ORT-SR1607	Ø6.0x70 mm
ORT-SR1608	Ø6.0x80 mm
ORT-SR1609	Ø6.0x90 mm
ORT-SR1610	Ø6.0x100 mm
ORT-SR1611	Ø6.0x110 mm
ORT-SR1612	Ø6.0x120 mm
ORT-SR1613	Ø6.0x130 mm
ORT-SR1614	Ø6.0x140 mm
ORT-SR1615	Ø6.0x150 mm
ORT-SR1616	Ø6.0x160 mm
ORT-SR1617	Ø6.0x170 mm
ORT-SR1618	Ø6.0x180 mm
ORT-SR1619	Ø6.0x190 mm
ORT-SR1620	Ø6.0x200 mm
ORT-SR1621	Ø6.0x210 mm
ORT-SR1622	Ø6.0x220 mm
ORT-SR1623	Ø6.0x230 mm
ORT-SR1624	Ø6.0x240 mm
ORT-SR1625	Ø6.0x250 mm
ORT-SR1626	Ø6.0x260 mm
ORT-SR1627	Ø6.0x270 mm
ORT-SR1628	Ø6.0x280 mm
ORT-SR1629	Ø6.0x290 mm
ORT-SR1630	Ø6.0x300 mm
ORT-SR1631	Ø6.0x310 mm
ORT-SR1632	Ø6.0x320 mm
ORT-SR1636	Ø6.0x360 mm
ORT-SR1640	Ø6.0x400 mm
ORT-SR1644	Ø6.0x440 mm
ORT-SR1648	Ø6.0x480 mm
ORT-SR1650	Ø6.0x500 mm
ORT-SR1660	Ø6.0x600 mm
ORT-SR3506	(3.5-6.0)x500 mm




POLAR

Spinal System 5.5/6.0 Rod

Colour Parameters as per Ø

						
Pink Ø 4.5	Purple Ø 5.0	Blue Ø 5.5	Turquoise Ø 6.0	Golden Ø 6.5	Ice Blue Ø 7.0	Green Ø 7.5



Features

- The Polar 5.5/6.0 Spinal System is low profile system.
- The screws feature a dual - lead thread design with self - tapping cuts, expediting advancement and reducing surgeon fatigue.
- Screws available in 4.5mm. 5.5mm. 6.5mm. 7.5mm and 8.5mm diameters and multiple lengths.
- Dual - lead thread pattern for faster insertation and increased pull out strength.
- Polyaxial Range of the motion for ease of use intraoperatively.
- Ability to accept diameter 5.5 and 6.0 mm rod.

Code	Size
ORT-KYA3520	3.5x20 mm
ORT-KYA3525	3.5x25 mm
ORT-KYA3530	3.5x30 mm
ORT-KYA3535	3.5x35 mm
ORT-KYA3540	3.5x40 mm
ORT-KYA3545	3.5x45 mm
ORT-KYA3550	3.5x50 mm
ORT-KYA3555	3.5x55 mm
ORT-KYA3560	3.5x60 mm
ORT-KYA4020	4.0x20 mm
ORT-KYA4025	4.0x25 mm
ORT-KYA4030	4.0x30 mm
ORT-KYA4035	4.0x35 mm
ORT-KYA4040	4.0x40 mm
ORT-KYA4045	4.0x45 mm
ORT-KYA4050	4.0x50 mm
ORT-KYA4055	4.0x55 mm
ORT-KYA4060	4.0x60 mm
ORT-KYA4520	4.5x20 mm
ORT-KYA4525	4.5x25 mm
ORT-KYA4530	4.5x30 mm
ORT-KYA4535	4.5x35 mm
ORT-KYA4540	4.5x40 mm
ORT-KYA4545	4.5x45 mm
ORT-KYA4550	4.5x50 mm
ORT-KYA4555	4.5x55 mm
ORT-KYA4560	4.5x60 mm
ORT-KYA5020	5.0x20 mm
ORT-KYA5025	5.0x25 mm
ORT-KYA5030	5.0x30 mm
ORT-KYA5035	5.0x35 mm
ORT-KYA5040	5.0x40 mm
ORT-KYA5045	5.0x45 mm
ORT-KYA5050	5.0x50 mm
ORT-KYA5055	5.0x55 mm
ORT-KYA5060	5.0x60 mm

Code	Size
ORT-KYA5520	5.5x20 mm
ORT-KYA5525	5.5x25 mm
ORT-KYA5530	5.5x30 mm
ORT-KYA5535	5.5x35 mm
ORT-KYA5540	5.5x40 mm
ORT-KYA5545	5.5x45 mm
ORT-KYA5550	5.5x50 mm
ORT-KYA5555	5.5x55 mm
ORT-KYA5560	5.5x60 mm
ORT-KYA6020	6.0x20 mm
ORT-KYA6025	6.0x25 mm
ORT-KYA6030	6.0x30 mm
ORT-KYA6035	6.0x35 mm
ORT-KYA6040	6.0x40 mm
ORT-KYA6045	6.0x45 mm
ORT-KYA6050	6.0x50 mm
ORT-KYA6055	6.0x55 mm
ORT-KYA6060	6.0x60 mm
ORT-KYA6520	6.5x20 mm
ORT-KYA6525	6.5x25 mm
ORT-KYA6530	6.5x30 mm
ORT-KYA6535	6.5x35 mm
ORT-KYA6540	6.5x40 mm
ORT-KYA6545	6.5x45 mm
ORT-KYA6550	6.5x50 mm
ORT-KYA6555	6.5x55 mm
ORT-KYA6560	6.5x60 mm
ORT-KYA7020	7.0x20 mm
ORT-KYA7025	7.0x25 mm
ORT-KYA7030	7.0x30 mm
ORT-KYA7035	7.0x35 mm
ORT-KYA7040	7.0x40 mm
ORT-KYA7045	7.0x45 mm
ORT-KYA7050	7.0x50 mm
ORT-KYA7055	7.0x55 mm
ORT-KYA7060	7.0x60 mm

Code	Size
ORT-KYA7520	7.5x20 mm
ORT-KYA7525	7.5x25 mm
ORT-KYA7530	7.5x30 mm
ORT-KYA7535	7.5x35 mm
ORT-KYA7540	7.5x40 mm
ORT-KYA7545	7.5x45 mm
ORT-KYA7550	7.5x50 mm
ORT-KYA7555	7.5x55 mm
ORT-KYA7560	7.5x60 mm
ORT-KYA8020	8.0x20 mm
ORT-KYA8025	8.0x25 mm
ORT-KYA8030	8.0x30 mm
ORT-KYA8035	8.0x35 mm
ORT-KYA8040	8.0x40 mm
ORT-KYA8045	8.0x45 mm
ORT-KYA8050	8.0x50 mm
ORT-KYA8055	8.0x55 mm
ORT-KYA8060	8.0x60 mm
ORT-KYA8065	8.0x65 mm
ORT-KYA8070	8.0x70 mm
ORT-KYA8075	8.0x75 mm
ORT-KYA8080	8.0x80 mm
ORT-KYA8085	8.0x85 mm
ORT-KYA8090	8.0x90 mm
ORT-KYA8095	8.0x95 mm
ORT-KYA80100	8.0x100 mm

Code	Size
ORT-KYA8520	8.5x20 mm
ORT-KYA8525	8.5x25 mm
ORT-KYA8530	8.5x30 mm
ORT-KYA8535	8.5x35 mm
ORT-KYA8540	8.5x40 mm
ORT-KYA8545	8.5x45 mm
ORT-KYA8550	8.5x50 mm
ORT-KYA8555	8.5x55 mm
ORT-KYA8560	8.5x60 mm
ORT-KYA8565	8.5x65 mm
ORT-KYA8570	8.5x70 mm
ORT-KYA8575	8.5x75 mm
ORT-KYA8580	8.5x80 mm
ORT-KYA8585	8.5x85 mm
ORT-KYA8590	8.5x90 mm
ORT-KYA8595	8.5x95 mm
ORT-KYA85100	8.5x100 mm



Code	Size
ORT-PPPS3525	Ø3.5x25 mm
ORT-PPPS3530	Ø3.5x30 mm
ORT-PPPS3535	Ø3.5x35 mm
ORT-PPPS3540	Ø3.5x40 mm
ORT-PPPS3545	Ø3.5x45 mm
ORT-PPPS4520	Ø4.5x20 mm
ORT-PPPS4525	Ø4.5x25 mm
ORT-PPPS4530	Ø4.5x30 mm
ORT-PPPS4535	Ø4.5x35 mm
ORT-PPPS4540	Ø4.5x40 mm
ORT-PPPS4545	Ø4.5x45 mm
ORT-PPPS5030	Ø5.0x30 mm
ORT-PPPS5035	Ø5.0x35 mm
ORT-PPPS5040	Ø5.0x40 mm
ORT-PPPS5045	Ø5.0x45 mm
ORT-PPPS5050	Ø5.0x50 mm
ORT-PPPS5520	Ø5.5x20 mm
ORT-PPPS5525	Ø5.5x25 mm
ORT-PPPS5530	Ø5.5x30 mm
ORT-PPPS5535	Ø5.5x35 mm
ORT-PPPS5540	Ø5.5x40 mm
ORT-PPPS5545	Ø5.5x45 mm
ORT-PPPS5550	Ø5.5x50 mm
ORT-PPPS5555	Ø5.5x55 mm
ORT-PPPS6035	Ø6.0x35 mm
ORT-PPPS6040	Ø6.0x40 mm
ORT-PPPS6045	Ø6.0x45 mm
ORT-PPPS6050	Ø6.0x50 mm
ORT-PPPS6055	Ø6.0x55 mm

Code	Size
ORT-PPPS6530	Ø6.5x30 mm
ORT-PPPS6535	Ø6.5x35 mm
ORT-PPPS6540	Ø6.5x40 mm
ORT-PPPS6545	Ø6.5x45 mm
ORT-PPPS6550	Ø6.5x50 mm
ORT-PPPS6555	Ø6.5x55 mm
ORT-PPPS7035	Ø7.0x35 mm
ORT-PPPS7040	Ø7.0x40 mm
ORT-PPPS7045	Ø7.0x45 mm
ORT-PPPS7050	Ø7.0x50 mm
ORT-PPPS7055	Ø7.0x55 mm
ORT-PPPS7530	Ø7.5x30 mm
ORT-PPPS7535	Ø7.5x35 mm
ORT-PPPS7540	Ø7.5x40 mm
ORT-PPPS7545	Ø7.5x45 mm
ORT-PPPS7550	Ø7.5x50 mm
ORT-PPPS7555	Ø7.5x55 mm
ORT-PPPS8030	Ø8.0x30 mm
ORT-PPPS8035	Ø8.0x35 mm
ORT-PPPS8040	Ø8.0x40 mm
ORT-PPPS8045	Ø8.0x45 mm
ORT-PPPS8050	Ø8.0x50 mm
ORT-PPPS8055	Ø8.0x55 mm
ORT-PPPS8060	Ø8.0x60 mm
ORT-PPPS8070	Ø8.0x70 mm
ORT-PPPS8080	Ø8.0x80 mm
ORT-PPPS8090	Ø8.0x90 mm
ORT-PPPS80100	Ø8.0x100 mm

Code	Size
ORT-PPPS8530	Ø8.5x30 mm
ORT-PPPS8535	8.5x35 mm
ORT-PPPS8540	8.5x40 mm
ORT-PPPS8545	8.5x45 mm
ORT-PPPS8550	8.5x50 mm
ORT-PPPS8555	8.5x55 mm
ORT-PPPS8560	8.5x60 mm
ORT-PPPS8570	8.5x70 mm
ORT-PPPS8580	8.5x80 mm
ORT-PPPS8590	8.5x90 mm
ORT-PPPS85100	8.5x100 mm





Code	Size
ORT-PPMAS3525	Ø3.5x25 mm
ORT-PPMAS3530	Ø3.5x30 mm
ORT-PPMAS3535	Ø3.5x35 mm
ORT-PPMAS3540	Ø3.5x40 mm
ORT-PPMAS3545	Ø3.5x45 mm
ORT-PPMAS4025	Ø4.0x25 mm
ORT-PPMAS4030	Ø4.0x30 mm
ORT-PPMAS4035	Ø4.0x35 mm
ORT-PPMAS4040	Ø4.0x40 mm
ORT-PPMAS4525	Ø4.5x25 mm
ORT-PPMAS4530	Ø4.5x30 mm
ORT-PPMAS4535	Ø4.5x35 mm
ORT-PPMAS4540	Ø4.5x40 mm
ORT-PPMAS4545	Ø4.5x45 mm
ORT-PPMAS5030	Ø5.0x30 mm
ORT-PPMAS5035	Ø5.0x35 mm
ORT-PPMAS5040	Ø5.0x40 mm
ORT-PPMAS5045	Ø5.0x45 mm
ORT-PPMAS5050	Ø5.0x50 mm
ORT-PPMAS5530	Ø5.5x30 mm
ORT-PPMAS5535	Ø5.5x35 mm
ORT-PPMAS5540	Ø5.5x40 mm
ORT-PPMAS5545	Ø5.5x45 mm
ORT-PPMAS5550	Ø5.5x50 mm
ORT-PPMAS5555	Ø5.5x55 mm

Code	Size
ORT-PPMAS6035	Ø6.0x35 mm
ORT-PPMAS6040	Ø6.0x40 mm
ORT-PPMAS6045	Ø6.0x45 mm
ORT-PPMAS6050	Ø6.0x50 mm
ORT-PPMAS6055	Ø6.0x55 mm
ORT-PPMAS6530	Ø6.5x30 mm
ORT-PPMAS6535	Ø6.5x35 mm
ORT-PPMAS6540	Ø6.5x40 mm
ORT-PPMAS6545	Ø6.5x45 mm
ORT-PPMAS6550	Ø6.5x50 mm
ORT-PPMAS6555	Ø6.5x55 mm
ORT-PPMAS7035	Ø7.0x35 mm
ORT-PPMAS7040	Ø7.0x40 mm
ORT-PPMAS7045	Ø7.0x45 mm
ORT-PPMAS7050	Ø7.0x50 mm
ORT-PPMAS7055	Ø7.0x55 mm
ORT-PPMAS7535	Ø7.5x35 mm
ORT-PPMAS7540	Ø7.5x40 mm
ORT-PPMAS7545	Ø7.5x45 mm
ORT-PPMAS7550	Ø7.5x50 mm
ORT-PPMAS7555	Ø7.5x55 mm
ORT-PPMAS8035	Ø8.0x35 mm
ORT-PPMAS8040	Ø8.0x40 mm
ORT-PPMAS8045	Ø8.0x45 mm
ORT-PPMAS8050	Ø8.0x50 mm
ORT-PPMAS8055	Ø8.0x55 mm

Polar Dual-Lead Spondylolisthesis Polyaxial Screw

Code	Size
ORT-PPRPS4530	Ø4.5x30 mm
ORT-PPRPS4535	Ø4.5x35 mm
ORT-PPRPS4540	Ø4.5x40 mm
ORT-PPRPS4545	Ø4.5x45 mm
ORT-PPRPS5530	Ø5.5x30 mm
ORT-PPRPS5535	Ø5.5x35 mm
ORT-PPRPS5540	Ø5.5x40 mm
ORT-PPRPS5545	Ø5.5x45 mm
ORT-PPRPS5550	Ø5.5x50 mm
ORT-PPRPS5555	Ø5.5x55 mm
ORT-PPRPS6035	Ø6.0x35 mm
ORT-PPRPS6040	Ø6.0x40 mm
ORT-PPRPS6045	Ø6.0x45 mm
ORT-PPRPS6050	Ø6.0x50 mm
ORT-PPRPS6535	Ø6.5x35 mm
ORT-PPRPS6540	Ø6.5x40 mm
ORT-PPRPS6545	Ø6.5x45 mm
ORT-PPRPS6550	Ø6.5x50 mm
ORT-PPRPS6555	Ø6.5x55 mm
ORT-PPRPS7035	Ø7.0x35 mm
ORT-PPRPS7040	Ø7.0x40 mm
ORT-PPRPS7045	Ø7.0x45 mm
ORT-PPRPS7050	Ø7.0x50 mm
ORT-PPRPS7055	Ø7.0x55 mm
ORT-PPRPS7535	Ø7.5x35 mm
ORT-PPRPS7540	Ø7.5x40 mm
ORT-PPRPS7545	Ø7.5x45 mm
ORT-PPRPS7550	Ø7.5x50 mm
ORT-PPRPS7555	Ø7.5x55 mm


Polar Dual-Lead Spondylolisthesis Monoaxial Screw

Code	Size
ORT-PPMRS4530	Ø4.5x30 mm
ORT-PPMRS4535	Ø4.5x35 mm
ORT-PPMRS4540	Ø4.5x40 mm
ORT-PPMRS4545	Ø4.5x45 mm
ORT-PPMRS5530	Ø5.5x30 mm
ORT-PPMRS5535	Ø5.5x35 mm
ORT-PPMRS5540	Ø5.5x40 mm
ORT-PPMRS5545	Ø5.5x45 mm
ORT-PPMRS5550	Ø5.5x50 mm
ORT-PPMRS6035	Ø6.0x35 mm
ORT-PPMRS6040	Ø6.0x40 mm
ORT-PPMRS6045	Ø6.0x45 mm
ORT-PPMRS6050	Ø6.0x50 mm
ORT-PPMRS6055	Ø6.0x55 mm
ORT-PPMRS6535	Ø6.5x35 mm
ORT-PPMRS6540	Ø6.5x40 mm
ORT-PPMRS6545	Ø6.5x45 mm
ORT-PPMRS6550	Ø6.5x50 mm
ORT-PPMRS6555	Ø6.5x55 mm
ORT-PPMRS7035	Ø7.0x35 mm
ORT-PPMRS7040	Ø7.0x40 mm
ORT-PPMRS7045	Ø7.0x45 mm
ORT-PPMRS7050	Ø7.0x50 mm
ORT-PPMRS7055	Ø7.0x55 mm
ORT-PPMRS7535	Ø7.5x35 mm
ORT-PPMRS7540	Ø7.5x40 mm
ORT-PPMRS7545	Ø7.5x45 mm
ORT-PPMRS7550	Ø7.5x50 mm
ORT-PPMRS7555	Ø7.5x55 mm



Polar Cannulated And Cemented Polyaxial Screw

Code	Size
ORT-PPCPS4530	Ø4.5x30 mm
ORT-PPCPS4535	Ø4.5x35 mm
ORT-PPCPS4540	Ø4.5x40 mm
ORT-PPCPS4545	Ø4.5x45 mm
ORT-PPCPS5530	Ø5.5x30 mm
ORT-PPCPS5535	Ø5.5x35 mm
ORT-PPCPS5540	Ø5.5x40 mm
ORT-PPCPS5545	Ø5.5x45 mm
ORT-PPCPS5550	Ø5.5x50 mm
ORT-PPCPS5555	Ø5.5x55 mm
ORT-PPCPS6530	Ø6.5x30 mm
ORT-PPCPS6535	Ø6.5x35 mm
ORT-PPCPS6540	Ø6.5x40 mm
ORT-PPCPS6545	Ø6.5x45 mm
ORT-PPCPS6550	Ø6.5x50 mm
ORT-PPCPS6555	Ø6.5x55 mm
ORT-PPCPS7530	Ø7.5x30 mm
ORT-PPCPS7535	Ø7.5x35 mm
ORT-PPCPS7540	Ø7.5x40 mm
ORT-PPCPS7545	Ø7.5x45 mm
ORT-PPCPS7550	Ø7.5x50 mm
ORT-PPCPS7555	Ø7.5x55 mm
ORT-PPCPS8035	Ø8.0x35 mm
ORT-PPCPS8040	Ø8.0x40 mm
ORT-PPCPS8045	Ø8.0x45 mm
ORT-PPCPS8050	Ø8.0x50 mm
ORT-PPCPS8055	Ø8.0x55 mm
ORT-PPCPS8060	Ø8.0x60 mm
ORT-PPCPS8070	Ø8.0x70 mm
ORT-PPCPS8080	Ø8.0x80 mm
ORT-PPCPS8090	Ø8.0x90 mm
ORT-PPCPS80100	Ø8.0x100 mm

Polar Cannulated And Cemented Monoaxial Screw

Code	Size
ORT-PPCMS5530	Ø5.5x30 mm
ORT-PPCMS5535	Ø5.5x35 mm
ORT-PPCMS5540	Ø5.5x40 mm
ORT-PPCMS5545	Ø5.5x45 mm
ORT-PPCMS5550	Ø5.5x50 mm
ORT-PPCMS5555	Ø5.5x55 mm
ORT-PPCMS6530	Ø6.5x30 mm
ORT-PPCMS6535	Ø6.5x35 mm
ORT-PPCMS6540	Ø6.5x40 mm
ORT-PPCMS6545	Ø6.5x45 mm
ORT-PPCMS6550	Ø6.5x50 mm
ORT-PPCMS6555	Ø6.5x55 mm
ORT-PPCMS7530	Ø7.5x30 mm
ORT-PPCMS7535	Ø7.5x35 mm
ORT-PPCMS7540	Ø7.5x40 mm
ORT-PPCMS7545	Ø7.5x45 mm
ORT-PPCMS7550	Ø7.5x50 mm
ORT-PPCMS7555	Ø7.5x55 mm



Polar Quad Lead Polyaxial Screw

Code	Size
ORT-PPMFS3525	Ø3.5x25 mm
ORT-PPMFS3530	Ø3.5x30 mm
ORT-PPMFS3535	Ø3.5x35 mm
ORT-PPMFS3540	Ø3.5x40 mm
ORT-PPMFS3545	Ø3.5x45 mm
ORT-PPMFS4520	Ø4.5x20 mm
ORT-PPMFS4525	Ø4.5x25 mm
ORT-PPMFS4530	Ø4.5x30 mm
ORT-PPMFS4535	Ø4.5x35 mm
ORT-PPMFS4540	Ø4.5x40 mm
ORT-PPMFS4545	Ø4.5x45 mm
ORT-PPMFS5030	Ø5.0x30 mm
ORT-PPMFS5035	Ø5.0x35 mm
ORT-PPMFS5040	Ø5.0x40 mm
ORT-PPMFS5045	Ø5.0x45 mm
ORT-PPMFS5050	Ø5.0x50 mm
ORT-PPMFS5520	Ø5.5x20 mm
ORT-PPMFS5525	Ø5.5x25 mm
ORT-PPMFS5530	Ø5.5x30 mm
ORT-PPMFS5535	Ø5.5x35 mm
ORT-PPMFS5540	Ø5.5x40 mm
ORT-PPMFS5545	Ø5.5x45 mm
ORT-PPMFS5550	Ø5.5x50 mm
ORT-PPMFS5555	Ø5.5x55 mm
ORT-PPMFS6035	Ø6.0x35 mm
ORT-PPMFS6040	Ø6.0x40 mm
ORT-PPMFS6045	Ø6.0x45 mm
ORT-PPMFS6050	Ø6.0x50 mm
ORT-PPMFS6055	Ø6.0x55 mm

Code	Size
ORT-PPMFS6530	Ø6.5x30 mm
ORT-PPMFS6535	Ø6.5x35 mm
ORT-PPMFS6540	Ø6.5x40 mm
ORT-PPMFS6545	Ø6.5x45 mm
ORT-PPMFS6550	Ø6.5x50 mm
ORT-PPMFS6555	Ø6.5x55 mm
ORT-PPMFS7035	Ø7.0x35 mm
ORT-PPMFS7040	Ø7.0x40 mm
ORT-PPMFS7045	Ø7.0x45 mm
ORT-PPMFS7050	Ø7.0x50 mm
ORT-PPMFS7055	Ø7.0x55 mm
ORT-PPMFS7530	Ø7.5x30 mm
ORT-PPMFS7535	Ø7.5x35 mm
ORT-PPMFS7540	Ø7.5x40 mm
ORT-PPMFS7545	Ø7.5x45 mm
ORT-PPMFS7550	Ø7.5x50 mm
ORT-PPMFS7555	Ø7.5x55 mm
ORT-PPMFS8030	Ø8.0x30 mm
ORT-PPMFS8035	Ø8.0x35 mm
ORT-PPMFS8040	Ø8.0x40 mm
ORT-PPMFS8045	Ø8.0x45 mm
ORT-PPMFS8050	Ø8.0x50 mm
ORT-PPMFS8055	Ø8.0x55 mm
ORT-PPMFS8070	Ø8.0x70 mm
ORT-PPMFS8080	Ø8.0x80 mm
ORT-PPMFS8090	Ø8.0x90 mm
ORT-PPMFS80100	Ø8.0x100 mm

Polar Quad Lead Spondylolisthesis Polyaxial Screw

Code	Size
ORT-PPMFRS5530	Ø5.5x30 mm
ORT-PPMFRS5535	Ø5.5x35 mm
ORT-PPMFRS5540	Ø5.5x40 mm
ORT-PPMFRS5545	Ø5.5x45 mm
ORT-PPMFRS5550	Ø5.5x50 mm
ORT-PPMFRS5555	Ø5.5x55 mm
ORT-PPMFRS6035	Ø6.0x35 mm
ORT-PPMFRS6040	Ø6.0x40 mm
ORT-PPMFRS6045	Ø6.0x45 mm
ORT-PPMFRS6050	Ø6.0x50 mm
ORT-PPMFRS6535	Ø6.5x35 mm
ORT-PPMFRS6540	Ø6.5x40 mm
ORT-PPMFRS6545	Ø6.5x45 mm
ORT-PPMFRS6550	Ø6.5x50 mm
ORT-PPMFRS6555	Ø6.5x55 mm
ORT-PPMFRS7035	Ø7.0x35 mm
ORT-PPMFRS7040	Ø7.0x40 mm
ORT-PPMFRS7045	Ø7.0x45 mm
ORT-PPMFRS7050	Ø7.0x50 mm
ORT-PPMFRS7055	Ø7.0x55 mm
ORT-PPMFRS7535	Ø7.5x35 mm
ORT-PPMFRS7540	Ø7.5x40 mm
ORT-PPMFRS7545	Ø7.5x45 mm
ORT-PPMFRS7550	Ø7.5x50 mm
ORT-PPMFRS7555	Ø7.5x55 mm

Polar Quad Lead Cannulated and Cemented Polyaxial Screw

Code	Size
ORT-PPMFC54530	Ø4.5x30 mm
ORT-PPMFC54535	Ø4.5x35 mm
ORT-PPMFC54540	Ø4.5x40 mm
ORT-PPMFC54545	Ø4.5x45 mm
ORT-PPMFC55530	Ø5.5x30 mm
ORT-PPMFC55535	Ø5.5x35 mm
ORT-PPMFC55540	Ø5.5x40 mm
ORT-PPMFC55545	Ø5.5x45 mm
ORT-PPMFC55550	Ø5.5x50 mm
ORT-PPMFC55555	Ø5.5x55 mm
ORT-PPMFC6035	Ø6.0x35 mm
ORT-PPMFC6040	Ø6.0x40 mm
ORT-PPMFC6045	Ø6.0x45 mm
ORT-PPMFC6050	Ø6.0x50 mm
ORT-PPMFC6535	Ø6.5x35 mm
ORT-PPMFC6540	Ø6.5x40 mm
ORT-PPMFC6545	Ø6.5x45 mm
ORT-PPMFC6550	Ø6.5x50 mm
ORT-PPMFC6555	Ø6.5x55 mm
ORT-PPMFC7035	Ø7.0x35 mm
ORT-PPMFC7040	Ø7.0x40 mm
ORT-PPMFC7045	Ø7.0x45 mm
ORT-PPMFC7050	Ø7.0x50 mm
ORT-PPMFC7055	Ø7.0x55 mm
ORT-PPMFC7535	Ø7.5x35 mm
ORT-PPMFC7540	Ø7.5x40 mm
ORT-PPMFC7545	Ø7.5x45 mm
ORT-PPMFC7550	Ø7.5x50 mm
ORT-PPMFC7555	Ø7.5x55 mm


Polar Quad Lead Spondylolisthesis Cannulated and Cemented Polyaxial Screw

Code	Size
ORT-MFRS5530	Ø5.5X30 mm
ORT-MFRS5535	Ø5.5X35 mm
ORT-MFRS5540	Ø5.5X40 mm
ORT-MFRS5545	Ø5.5X45 mm
ORT-MFRS5550	Ø5.5X50 mm
ORT-MFRS5555	Ø5.5X55 mm
ORT-MFRS6035	Ø6.0X35 mm
ORT-MFRS6040	Ø6.0X40 mm
ORT-MFRS6045	Ø6.0X45 mm
ORT-MFRS6050	Ø6.0X50 mm
ORT-MFRS6535	Ø6.5X35 mm
ORT-MFRS6540	Ø6.5X40 mm
ORT-MFRS6545	Ø6.5X45 mm
ORT-MFRS6550	Ø6.5X50 mm
ORT-MFRS6555	Ø6.5X55 mm
ORT-MFRS7035	Ø7.0X35 mm
ORT-MFRS7040	Ø7.0X40 mm
ORT-MFRS7045	Ø7.0X45 mm
ORT-MFRS7050	Ø7.0X50 mm
ORT-MFRS7055	Ø7.0X55 mm
ORT-MFRS7535	Ø7.5X35 mm
ORT-MFRS7540	Ø7.5X40 mm
ORT-MFRS7545	Ø7.5X45 mm
ORT-MFRS7550	Ø7.5X50 mm
ORT-MFRS7555	Ø7.5X55 mm



M.I.S

Code	Size
ORT-SPA3525	Ø3.5x25 mm
ORT-SPA3530	Ø3.5x30 mm
ORT-SPA3535	Ø3.5x35 mm
ORT-SPA3540	Ø3.5x40 mm
ORT-SPA3545	Ø3.5x45 mm
ORT-SPA4520	Ø4.5x20 mm
ORT-SPA4525	Ø4.5x25 mm
ORT-MFS5030	Ø5.0x30 mm
ORT-MFS5035	Ø5.0x35 mm
ORT-MFS5040	Ø5.0x40 mm
ORT-MFS5045	Ø5.0x45 mm
ORT-MFS5050	Ø5.0x50 mm
ORT-SPA5520	Ø5.5x20 mm
ORT-MFS5525	Ø5.5x25 mm
ORT-MFS5530	Ø5.5x30 mm
ORT-MFS5535	Ø5.5x35 mm
ORT-MFS5540	Ø5.5x40 mm
ORT-MFS5545	Ø5.5x45 mm
ORT-MFS5550	Ø5.5x50 mm
ORT-MFS5555	Ø5.5x55 mm



Code	Size
ORT-MFS6035	Ø6.0x35 mm
ORT-MFS6040	Ø6.0x40 mm
ORT-MFS6045	Ø6.0x45 mm
ORT-MFS6050	Ø6.0x50 mm
ORT-MFS6055	Ø6.0x55 mm
ORT-MFS6530	Ø6.5x30 mm
ORT-MFS6535	Ø6.5x35 mm
ORT-MFS6540	Ø6.5x40 mm
ORT-MFS6545	Ø6.5x45 mm
ORT-MFS6550	Ø6.5x50 mm
ORT-MFS6555	Ø6.5x55 mm
ORT-MFS7035	Ø7.0x35 mm
ORT-MFS7040	Ø7.0x40 mm
ORT-MFS7045	Ø7.0x45 mm
ORT-MFS7050	Ø7.0x50 mm
ORT-MFS7055	Ø7.0x55 mm

Code	Size
ORT-MFS7530	Ø7.5x30 mm
ORT-MFS7535	Ø7.5x35 mm
ORT-MFS7540	Ø7.5x40 mm
ORT-MFS7545	Ø7.5x45 mm
ORT-MFS7550	Ø7.5x50 mm
ORT-MFS7555	Ø7.5x55 mm
ORT-MFS8030	Ø8.0x30 mm
ORT-MFS8035	Ø8.0x35 mm
ORT-MFS8040	Ø8.0x40 mm
ORT-MFS8045	Ø8.0x45 mm
ORT-MFS8050	Ø8.0x50 mm
ORT-MFS8055	Ø8.0x55 mm
ORT-MFS8070	Ø8.0x70 mm
ORT-MFS8080	Ø8.0x80 mm
ORT-MFS8090	Ø8.0x90 mm
ORT-MFS80100	Ø8.0x100 mm

Rod 5.5

Code	Size
ORT-SR1504	40 mm
ORT-SR1505	50 mm
ORT-SR1506	60 mm
ORT-SR1507	70 mm
ORT-SR1508	80 mm
ORT-SR1509	90 mm
ORT-SR1510	100 mm
ORT-SR1511	110 mm
ORT-SR1512	120 mm
ORT-SR1513	130 mm
ORT-SR1514	140 mm
ORT-SR1515	150 mm
ORT-SR1515	160 mm
ORT-SR1517	170 mm
ORT-SR1518	180 mm
ORT-SR1519	190 mm
ORT-SR1520	200 mm
ORT-SR1521	210 mm
ORT-SR1522	220 mm
ORT-SR1523	230 mm
ORT-SR1524	240 mm
ORT-SR1525	250 mm
ORT-SR1526	260 mm
ORT-SR1527	270 mm
ORT-SR1528	280 mm
ORT-SR1529	290 mm
ORT-SR1530	300 mm
ORT-SR1531	310 mm
ORT-SR1532	320 mm
ORT-SR1536	360 mm
ORT-SR1540	400 mm
ORT-SR1544	440 mm
ORT-SR1548	480 mm
ORT-SR1550	500 mm
ORT-SR1560	600 mm

Type 2 Rod 5.5

Code	Size
ORT-MR1504	40 mm
ORT-MR1505	50 mm
ORT-MR1506	60 mm
ORT-MR1507	70 mm
ORT-MR1508	80 mm
ORT-MR1509	90 mm
ORT-MR1510	100 mm
ORT-MR1511	110 mm
ORT-MR1512	120 mm
ORT-MR2504	40 mm
ORT-MR2505	50 mm
ORT-MR2506	60 mm
ORT-MR2507	70 mm
ORT-MR2508	80 mm
ORT-MR2509	90 mm
ORT-MR2510	100 mm
ORT-MR2511	110 mm
ORT-MR2512	120 mm

Straight

Curved

Transverse Connector

Code	Size
ORT-TL140	40 mm
ORT-TL150	50 mm
ORT-TL160	60 mm
ORT-TL170	70 mm
ORT-TL180	80 mm

Code	Size
ORT-TLR140	40 mm
ORT-TLR150	50 mm
ORT-TLR160	60 mm
ORT-TLR170	70 mm
ORT-TLR180	80 mm
ORT-TLH	-

Multiaxial Transverse Connector

Code	Size
ORT-MTL2040	20-40 mm
ORT-MTL4060	40-60 mm
ORT-MTL6080	60-80 mm



Lateral Connector

Code	Size
ORT-LCNT15	15 mm
ORT-LCNT20	20 mm
ORT-LCNT25	25 mm
ORT-LCNT30	30 mm



Lateral Open Connector

Code	Size
ORT-LONT15	15 mm
ORT-LONT20	20 mm
ORT-LONT25	25 mm
ORT-LONT30	30 mm



Pedicular Hooks

Code	Size
ORT-PH0805	8x5 mm
ORT-PH0807	8x7 mm
ORT-PH0809	8x9 mm
ORT-PHF0505	5x5 mm
ORT-PHF0507	5x7 mm



Set Screw

Code

ORT-SSS



Lateral Connector

Code	Size
ORT-LH0505	5x5 mm
ORT-LH0507	5x7 mm
ORT-LH0509	5x9 mm
ORT-LH0706	7x6 mm
ORT-LH0707	7x7 mm
ORT-LH0709	7x9 mm
ORT-LH0711	7x11 mm
ORT-LHF0505	5x5 mm
ORT-LHF0507	5x7 mm
ORT-LHLA709	7x9 mm
ORT-LHLA711	7x11 mm
ORT-LHRA709	7x9 mm
ORT-LHRA711	7x11 mm



Extension Connector For 1 Rod

Code

ORT-ECNT01



Extension Connector For 2 Rod

Code

ORT-ECNT02



Open Extension Connector For 2 Rod Small

Code

ORT-ECNTAK



Open Extension Connector For 2 Rod Large

Code

ORT-ECNTAL



Lumbar Corpectomy Cage



Features

- Ti-6Al-4V ELI Titanium alloy.
- Anatomical design.
- Efficient design providing wide grafting space.
- Controlled in-situ expansion.
- Easy instrumentation.






Code	Dia	x	Closed Lenght	x	Open Lenght
ORTVTCD202534	20	x	25	x	34
ORTVTCD203244	20	x	32	x	44
ORTVTCD204258	20	x	42	x	58
ORTVTCD205684	20	x	56	x	84


Code	Dia	x	Closed Lenght	x	Open Lenght
ORTVTCD242329	24	x	23	x	29
ORTVTCD242840	24	x	28	x	40
ORTVTCD243854	24	x	38	x	54
ORTVTCD245276	24	x	52	x	76
ORTVTCD2474114	24	x	74	x	114

Code	Dia	x	Closed Lenght	x	Open Lenght
ORTVTCD283044	28	x	30	x	44
ORTVTCD284258	28	x	42	x	58
ORTVTCD285684	28	x	56	x	84
ORTVTCD2882130	28	x	82	x	130
ORTVTCD28100140	28	x	100	x	140



Code	Dia	x	Closed Lenght	x	Open Lenght	Angled
ORTVTCA202534	20	x	25	x	34	6°
ORTVTCA203244	20	x	32	x	44	6°
ORTVTCA204258	20	x	42	x	58	6°
ORTVTCA205684	20	x	56	x	84	6°

Code	Dia	x	Closed Lenght	x	Open Lenght	Angled
ORTVTCA242329	24	x	23	x	29	6°
ORTVTCA242840	24	x	28	x	40	6°
ORTVTCA243854	24	x	38	x	54	6°
ORTVTCA245276	24	x	52	x	76	6°
ORTVTCA2474114	24	x	74	x	114	6°




Code	Dia	x	Closed Lenght	x	Open Lenght	Angled
ORTVTCA283044	28	x	30	x	44	6°
ORTVTCA284258	28	x	42	x	58	6°
ORTVTCA285684	28	x	56	x	84	6°
ORTVTCA2882130	28	x	82	x	130	6°
ORTVTCA28100140	28	x	100	x	140	6°



- **Simplicity of use.**
- **Dynamic fixation and stabilisation.**
- **Reduces stresses at the facet joints.**
- **Opens up the foramina.**
- **Prevent upfolding ligamentum flavum.**

Code	Size
ORT-8UDPL	08 mm
ORT-10UDPL	10 mm
ORT-12UDPL	12 mm
ORT-14UDPL	14 mm
ORT-16UDPL	16 mm


Tlif Peek Cage



Features

- Bi-convex graft contour, with a 4° angle of lordosis, conforms easily to the concavity of the endplate.
- Structure on the surface, minimizing the risk of expulsion.
- Large contact area, minimizing the risk of subsidence.
- Large graft area, maximizing the chances of a successful fusion.
- Simple, ergonomic and intelligent instrumentation.





Code	Height	Lenght	Width
ORT-BC2407	7	24	10
ORT-BC2408	8	24	10
ORT-BC2409	9	24	10
ORT-BC2410	10	24	10
ORT-BC2411	11	24	10
ORT-BC2412	12	24	10
ORT-BC2413	13	24	10

Code	Height	Lenght	Width	Angled
ORT-BCA2407	7	24	10	4°
ORT-BCA2408	8	24	10	4°
ORT-BCA2409	9	24	10	4°
ORT-BCA2410	10	24	10	4°
ORT-BCA2411	11	24	10	4°
ORT-BCA2412	12	24	10	4°
ORT-BCA2413	13	24	10	4°

Code	Height	Lenght	Width	Angled
ORT-BCA2807	7	28	10	4°
ORT-BCA2808	8	28	10	4°
ORT-BCA2809	9	28	10	4°
ORT-BCA2810	10	28	10	4°
ORT-BCA2811	11	28	10	4°
ORT-BCA2812	12	28	10	4°
ORT-BCA2813	13	28	10	4°


Code	Height	Lenght	Width	Angled
ORT-BCA3207	7	32	10	4°
ORT-BCA3208	8	32	10	4°
ORT-BCA3209	9	32	10	4°
ORT-BCA3210	10	32	10	4°
ORT-BCA3211	11	32	10	4°
ORT-BCA3212	12	32	10	4°
ORT-BCA3213	13	32	10	4°

Code	Height	Lenght	Width
ORT-BC2807	7	28	10
ORT-BC2808	8	28	10
ORT-BC2809	9	28	10
ORT-BC2810	10	28	10
ORT-BC2811	11	28	10
ORT-BC2812	12	28	10
ORT-BC2813	13	28	10

Code	Height	Lenght	Width
ORT-BC3207	7	32	10
ORT-BC3208	8	32	10
ORT-BC3209	9	32	10
ORT-BC3210	10	32	10
ORT-BC3211	11	32	10
ORT-BC3212	12	32	10
ORT-BC3213	13	32	10




Plif Peek Cage



Features

- Self-distracting, bullet-shaped nose, facilitating insertion.
- Convex-shaped Implants, designed to the patient anatomy and to allow more accurate sizing.
- Teeth on the surface, minimizing the risk of expulsion.
- Large contact area, minimizing the risk of subsidence.
- Simple, ergonomic and intelligent instrumentation.





Code	Height	Lenght	Width	Angled
ORT-LCA2206	6	22	10	4°
ORT-LCA2207	7	22	10	4°
ORT-LCA2208	8	22	10	4°
ORT-LCA2209	9	22	10	4°
ORT-LCA2210	10	22	10	4°
ORT-LCA2211	11	22	10	4°
ORT-LCA2212	12	22	10	4°
ORT-LCA2214	14	22	10	4°
ORT-LCA2216	16	22	10	4°

Code	Height	Lenght	Width	Angled
ORT-LCA2406	6	24	10	4°
ORT-LCA2407	7	24	10	4°
ORT-LCA2408	8	24	10	4°
ORT-LCA2409	9	24	10	4°
ORT-LCA2410	10	24	10	4°
ORT-LCA2411	11	24	10	4°
ORT-LCA2412	12	24	10	4°




Code	Height	Lenght	Width
ORT-LCD2206	6	22	10
ORT-LCD2207	7	22	10
ORT-LCD2208	8	22	10
ORT-LCD2209	9	22	10
ORT-LCD2210	10	22	10
ORT-LCD2211	11	22	10
ORT-LCD2212	12	22	10

Code	Height	Lenght	Width
ORT-LCD2406	6	24	10
ORT-LCD2407	7	24	10
ORT-LCD2408	8	24	10
ORT-LCD2409	9	24	10
ORT-LCD2410	10	24	10
ORT-LCD2411	11	24	10
ORT-LCD2412	12	24	10



Expandable Plif Peek Cage



Features

- Self-distracting, bulleted nose, facilitating insertion.
- Easy to implant with sharp ended design.
- Teeth on the surface, minimizing the risk of expulsion.
- Graft area, maximizing the chances of a successful fusion.
- Large contact area, minimizing the risk of subsidence.





Code	Height	Lenght	Width
ORT-EPC1007	7	24	10
ORT-EPC1008	8	24	10
ORT-EPC1009	9	24	10
ORT-EPC1010	10	24	10
ORT-EPC1011	11	24	10
ORT-EPC1012	12	24	10
ORT-EPC1013	13	24	10






Pedis Spinal System

PEDIS

Spinal System 5.0mm Rod

Colour Parameters as per Ø

								
Light Yellow Ø 3.5	Pink Ø 4.0	Yellow Ø 4.5	Turquoise Ø 5.0	Light Blue Ø 5.5	Blue Ø 6.0	Purple Ø 6.5	Light Green Ø 7.0	Green Ø 7.5



Features

- The pedis pediatric pedicle system consists of polyaxial screws, and reduction screw, various rods, and locking set screws to provide efficient and secure top-loading, rigid fixation in posterior lumbar fusions for pediatric patient.

Pediatric Polyaxial Screw



Code	Size
ORT-PPAS3520	Ø3.5x20 mm
ORT-PPAS3525	Ø3.5x25 mm
ORT-PPAS3530	Ø3.5x30 mm
ORT-PPAS3535	Ø3.5x35 mm
ORT-PPAS3540	Ø3.5x40 mm
ORT-PPAS4020	Ø4.0x20 mm
ORT-PPAS4025	Ø4.0x25 mm
ORT-PPAS4030	Ø4.0x30 mm
ORT-PPAS4035	Ø4.0x35 mm
ORT-PPAS4040	Ø4.0x40 mm
ORT-PPAS4530	Ø4.5x30 mm
ORT-PPAS4535	Ø4.5x35 mm
ORT-PPAS4540	Ø4.5x40 mm
ORT-PPAS4545	Ø4.5x45 mm
ORT-PPAS5030	Ø5.0x30 mm
ORT-PPAS5035	Ø5.0x35 mm
ORT-PPAS5040	Ø5.0x40 mm
ORT-PPAS5045	Ø5.0x45 mm
ORT-PPAS5530	Ø5.5x30 mm
ORT-PPAS5535	Ø5.5x35 mm
ORT-PPAS5540	Ø5.5x40 mm
ORT-PPAS5545	Ø5.5x45 mm
ORT-PPAS5550	Ø5.5x50 mm
ORT-PPAS6030	Ø6.0x30 mm
ORT-PPAS6035	Ø6.0x35 mm
ORT-PPAS6040	Ø6.0x40 mm
ORT-PPAS6045	Ø6.0x45 mm
ORT-PPAS6050	Ø6.0x50 mm
ORT-PPAS6530	Ø6.5x30 mm
ORT-PPAS6535	Ø6.5x35 mm
ORT-PPAS6540	Ø6.5x40 mm
ORT-PPAS6545	Ø6.5x45 mm
ORT-PPAS6550	Ø6.5x50 mm
ORT-PPAS7030	Ø7.0x30 mm
ORT-PPAS7035	Ø7.0x35 mm
ORT-PPAS7040	Ø7.0x40 mm
ORT-PPAS7045	Ø7.0x45 mm
ORT-PPAS7050	Ø7.0x50 mm
ORT-PPAS7055	Ø7.0x55 mm
ORT-PPAS7060	Ø7.0x60 mm

Pediatric Polyaxial Spondylolisthesis Screw



Code	Size
ORT-SPPAS4035	Ø4.0x35 mm
ORT-SPPAS4040	Ø4.0x40 mm
ORT-SPPAS4045	Ø4.0x45 mm
ORT-SPPAS5035	Ø5.0x35 mm
ORT-SPPAS5040	Ø5.0x40 mm
ORT-SPPAS5045	Ø5.0x45 mm
ORT-SPPAS5535	Ø5.5x35 mm
ORT-SPPAS5540	Ø5.5x40 mm
ORT-SPPAS5545	Ø5.5x45 mm
ORT-SPPAS6530	Ø6.5x30 mm
ORT-SPPAS6535	Ø6.5x35 mm
ORT-SPPAS6540	Ø6.5x40 mm
ORT-SPPAS6545	Ø6.5x45 mm
ORT-SPPAS6550	Ø6.5x50 mm

Pediatric Laminar Hooks

Code	Size
ORT-PLH0455	4.0x5.5 mm
ORT-PLH0465	4.0x6.5 mm
ORT-PLH0475	4.0x7.5 mm
ORT-PLH0555	5.0x5.5 mm
ORT-PLH0565	5.0x6.5 mm
ORT-PLH0575	5.0x7.5 mm
ORT-PLH0607	6.0x7.0 mm



Pediatric Rod 5.0

Code	Size
ORT-PPR40	Ø5.0x40 mm
ORT-PPR60	Ø5.0x60 mm
ORT-PPR80	Ø5.0x80 mm
ORT-PPR100	Ø5.0x100 mm
ORT-PPR120	Ø5.0x120 mm
ORT-PPR140	Ø5.0x140 mm
ORT-PPR160	Ø5.0x160 mm
ORT-PPR180	Ø5.0x180 mm
ORT-PPR200	Ø5.0x200 mm
ORT-PPR250	Ø5.0x250 mm
ORT-PPR300	Ø5.0x300 mm
ORT-PPR400	Ø5.0x400 mm
ORT-PPR500	Ø5.0x500 mm



Pediatric Pedicular Hooks

Code	Size
ORT-PPH0605	6.0x5.0 mm
ORT-PPH0705	7.0x5.0 mm
ORT-PPH0707	7.0x7.0 mm
ORT-PPH7555	7.5x5.5 mm
ORT-PPH7565	7.5x6.5 mm
ORT-PPH0905	9.0x5.0 mm
ORT-PPH0907	9.0x7.0 mm



Pediatric Transverse Connector

Code	Size
ORT-PTC40	40 mm
ORT-PTC50	50 mm
ORT-PTC60	60 mm



Pediatric Domino Connector

Code	Size
ORT-PPAK	S
ORT-PPAO	M
ORT-PPAB	L



Pediatric Rod Connector

Code	Size
ORT-PPRC5060	5.0/6.0
ORT-PPRC7550	5.0/5.0



Pediatric Set Screw

Code	
ORT-PSS	

Pediatric Lateral Connector

Code	
ORT-PPLC	

Pediatric Extension Connector

Code	
ORT-OCNTPE	